# PCRFast<sup>®</sup> Haselnuss

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Lit.: Holzhauser T, Stephan O, Vieths S

„Detection of potentially allergenic hazelnut (*Corylus avellana*) residues in food: A comparative study with DNA PCR-ELISA and protein sandwich ELISA“, J Agric Food Chem. (2002) 50, 5808-5815

Lit: Sambrook J, Fritsch E F, Maniatis T

„Molecular cloning. A laboratory manual“, (2<sup>nd</sup> ed.), New York: Cold Spring Harbor Laboratory Press.

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09/07/2008

# PCRFast<sup>®</sup> Hazelnut

Realtime (SYBR<sup>®</sup> Green) and gel detection

## PCR - test for the detection of hazelnut

Type No.: IF / AG1001

In-vitro test

Storage: 2 - 8 °C (35.6 - 46.4 °F)

### Brief information

Easy-to-use molecular biological test (PCR) for the detection of hazelnut (*Corylus avellana*, detection of the "major allergen protein" gene, 78bp) in food and pharmaceutical products. The test kit contains 96 reaction vials. All reaction vials contain a specific primer pair, and in addition 48 reaction vials (red-coloured) contain specific hazelnut DNA for the PCR positive control and for the inhibition control.

### Usage

For the general usage of the kit see "PCRFast<sup>®</sup> - General Information". Please read the instruction carefully.

### Cycler profile

10 min	95 °C (203 °F)	<b>45 cycles</b>
15 sec	95 °C (203 °F)	
60 sec	62 °C (143.6 °F)	

The sequential identity of the amplified DNA (78bp) could be verified by a restriction analysis with *Bam*H I. This results in two fragments with a length of 20bp and 58bp.

### Sensitivity

The limit of detection is < 10 copies.

### Specificity

PCRFast<sup>®</sup> Hazelnut is 100 % specific for hazelnut. The following species have been tested with 100 ng DNA for cross reactivity:

species		species		species		species	
hazelnut	+	pecan	-	celery	-	barley	-
peanut	-	pistachio	-	cabbage	-	oat	-
almond	-	brazil nut	-	orange	-	spelt	-
cashew	-	apricot	-	tangerine	-	buckwheat	-
macadamia	-	soy	-	wheat	-		
walnut	-	corn	-	rye	-		

Tab. 1: Specificity PCRFast<sup>®</sup> Hazelnut

+ : band with 78bp

- : no band with 78bp

Lit.: Holzhauser T, Stephan O, Vieths S

"Detection of potentially allergenic hazelnut (*Corylus avellana*) residues in food: A comparative study with DNA PCR-ELISA and protein sandwich ELISA", J Agric Food Chem. (2002) 50, 5808-5815

Lit: Sambrook J, Fritsch E F, Maniatis T

„Molecular cloning. A laboratory manual“, (2<sup>nd</sup>ed.), New York: Cold Spring Harbor Laboratory Press.

### Certificate of Quality

PCRFast<sup>®</sup> Hazelnut,  
complies specification.

Lot: **THA 11111**  
quality assurance: A.Schuhr

09/07/2008