PCRFast[®] Weizen / Gerste / Roggen

Realtime (Sonde)

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PCRFast[®] Wheat / Barley / Rye

Realtime (probe)

PCR - test for the detection of wheat, barley and rye

Type No.: IF / AR1009

In-vitro test

Storage: 2 - 8 °C (35.6 - 46.4 °F)

Brief information

Easy-to-use molecular biological test (PCR) for the detection of wheat (*triticum aestivum*), barley (*Hordeum vulgare*) and rye (*Secale cereale*), probe with **Fam** reporter and non fluorescence quencher) in food, feeding stuff and pharmaceutical products (detection of the chloroplast trnL gene). The test kit has 96 reaction vials, each reaction vial contains a specific primer pair, the probe and an internal amplification control (probe with **HEX** reporter and non fluorescence quencher) for investigating possible inhibiting effects (ITC). 8 reaction vials (red marking) contain additionally cereal DNA for the PCR positive controls.

Usage

For the general usage of the kit see "PCRFast[®] - General Information". Please read the instruction carefully.

Cycler profile

10 min	95 °C (203 °F)	45 cycles
15 sec	95 °C (203 °F)	
60 sec	55 °C (131 °F)	
60 sec	72 °C (161.6 °F)	

Sensitivity

The limit of detection is < 10 copies.

Specificity

PCRFast[®] Wheat/Barley/Rye is 100 % specific for wheat, barley and rye. The following species have been tested with 100 ng DNA for cross reactivity:

species		species		species		species	
wheat	+	barley	+	rye	+	corn	-
rice	-	soy	-	sesame	-	kidney bean	-
sunflowers	-						

Tab. 1: Specificity PCRFast[®] Wheat/Barley/Rye

+ : amplification

- : no amplification

Lit.: Dahinden I, von Büren M, Lüthy J

„A quantitative competitive PCR system to detect contamination of wheat, barley or rye in gluten-free food for coeliac patients“ (2001) J. Food Res: Technol 212,228-233

Lit: Sambrook J, Fritsch E F, Maniatis T

„Molecular cloning. A laboratory manual“, (3rded.), New York: Cold Spring Harbor Laboratory Press.

Certificate of Quality

PCRFast[®] Wheat/Barley/Rye,
complies specification.

Lot: **RGE 1111**
quality assurance: A.Schuhr

18/07/2008