**PCRFast[®]
Norovirus**

**Realtime (Sonde)
Realtime (probe)**

IF/MR1022

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Kurzinformation

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Hoehne M., Schreier E. “Detection of Norovirus genogroup I and II by multiplex real-time RT-PCR using a 3'-minor groove binder-DNA probe” BMC Infectious Diseases (2006), **6 (69)**.

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PCRFast[®] Norovirus Realtime (probe)

Brief information

Simple, molecular biological test (PCR) for the qualitative detection of norovirus in realtime (probe with Fam reporter and non fluorescence quencher) in foodstuffs (detection of the genotype I and genotype II). The test kit has 96 reaction vials. Each reaction vial contain specific primer pairs and probes for the analysis of samples and negative controls (NTC). Additionally every reaction tube contains a reverse transcription control (RTC, probe with HEX reporter and non fluorescence quencher) for testing the reverse transcription and investigating possible inhibition effects.

8 reaction vials (red marking) contain additionally norovirus - cDNA for the positive controls (PTC).

Time required:	RNA - Extraction.....	45 min
(10 samples)	PCR setup.....	15 min
	PCR.....	2.0 h

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PCRFast[®]

is a registered trademark of the ifp Institut für Produktqualität GmbH.
ifp also offers contract analysis.

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1. Principle of the test

PCRFast[®] Norovirus is a simple, molecular biological test (PCR) for the qualitative detection of norovirus (genotype I and genotype II) in food. The analytical procedure described here complies with the international standards (e.g. ISO/DIS 20838) for PCR analysis.

All the reaction vials contain the optimum amount of the specific primers and probes required for the PCR reaction and a **reverse transcription control (RTC)** consisting of MS2 phage specific primers and probe (HEX labelled probe with non fluorescence quencher) as well as MS2 phage RNA for checking the reverse transcription and inhibition control.

The **88 colourless reaction vials (coulerless)** contain specific primer pairs and probes (FAM labelled with non fluorescence quencher) used for detecting **norovirus RNA by reverse transcription and subsequent amplification of the norovirus cDNA** in the extracted sample and for the **negative controls NTC** (checking the MasterMix for contamination).

The **8 red marked reaction vials** contain **norovirus- cDNA** as well as the specific primers and probes (Fam labelled with non fluorescence quencher). These are used for the **positive controls PTC** (checking the functionality of MasterMix).

Place **12.75 µl of MasterMix** into the reaction vial (12.5 µl 2x QuantiTect Probe RT-PCR Master Mix and 0.25 µl QuantiTect RT Mix per reaction) and add **12.25 µl of the RNA sample extract**. The MasterMix already contains polymerase, transcriptase, nucleotides and magnesium chloride in optimum concentrations. The target will be transcribed and amplified in the PCR thermocycler. Is the target present the probe hybridises on the synthesised DNA fragments. Thereby the reporter (Fam / HEX) and the quencher (non fluorescence) of the probe will be separated and the reporter radiates a fluorescence signal (Fam: 520 nm / HEX: 553 nm). The intensity of the signal increases with the amount of synthesised DNA which can be measured and displayed. If the target sequence is not present no increase in the fluorescence is detectable.

Note:

RNA solutions are very highly sensitive to ribonucleases. Concentrated working and the avoiding of standing times of RNA-solutions (RNA extracts, PCR-setup;...) are crucial for a sensitive detection of Norovirus. To avoid longer standing times it is possible to minimize the ribonuklease activity by freezing the PCR – setup before the amplification at no less than – 20°C (- 4 °F).

2. Package contents

The test kit contains:

- 11 x **strips (colourless)** with 8 x 0.2 ml reaction vials each (88 reaction vials), coated with specific primers and probes (for sample determinations and for the negative controls, NTC) and containing a reverse transcription control (RTC, checking the reverse transcription and inhibition effects).
- 1 x **strips (red)** with 8 x 0.2 ml reaction vials each (8 reaction vials), coated with specific primers, probes and additionally Norovirus - cDNA (for the positive controls PTC) and containing a reverse transcription control (RTC; checking the reverse transcription and inhibition effects).

3. Additionally required instruments and reagents

Lab material and instruments

- vortexer
- biocentrifuge, min. 14,000 x g (for 1.5 ml and 2.0 ml reaction vials)
- PCR thermocycler:**
 - e.g. Applied Biosystems 7500, STRATAGENE Mx 3005 P, Mx 3000 P
- microliter pipettes 2 - 20 µl, 20 - 200 µl, 100 - 1,000 µl, e.g. Gilson Pipetman P with filter tips
- 1.5 ml or 2.0 ml reaction vials
- gauze swab, plastic shaft recommended

Reagents

- ethanol >95 % (denatured)
- 0.1 x TE - buffer (TRIS 1 mmol / l, EDTA 0.1 mmol / l)
- water deionised, DNA free
- DNA-extraction kit, e.g. Qiagen „QIAamp[®] Viral RNA Kit

PCR

Realtime

- double concentrated MasterMix (reverse transcription and PCR)
QuantiTect™ Probe RT - PCR Kit, Qiagen, No.: 204443, Brilliant® II
QRT-PCR MasterMix Kit, Stratagene, No.:600809

4. Precautions

- all tasks should be performed taking all precautions common in labs, and all PCR tasks should be performed in accordance with the CEN / ISO recommendations
- to avoid carryovers, perform all work wearing protective gloves and use filter tips
- perform sample preparation, PCR setup and detection in separate rooms

5. Storage instructions

Store the test kit and the reaction vials at 2 - 8 °C (35.6 - 46.4 °F).

6. RNA Isolation

Food samples are extracted according to § 64 of the German Foodstuffs and Animal Feed Code (LFGB, L00.00 - 112).

Drop off the norovirus with a steril, moistened (e.g. PBS) swab from the food material. Transfer the viruses directly in the Lysesbuffer (in commercial RNA – Extraction kits included, see L00.00 - 112). Proceed in comply with manufacturer description (e.g. Qiagen “Qiamp® Viral RNA Kit”)

With each analyzing series one sample without matrix should be done as an extraction control. A positive extraction control indicates a contamination in the reagents.

7. PCR setup

7.1. Preparations

- remove strip from the plastic bag and separate required number of reaction vials
- place the remaining strips / reaction vials and the desiccant back into the bag, seal bag and store at 2 - 8 °C (35.6 - 46.4 °F)
- have the required amount of MasterMix at hand (12.75 µl per reaction)

7.2. Procedure

–pipette the following volumes into the reaction vials:

	reaction vial	Master-Mix	extracted sample RNA
sample	colourless	12.75 µl	12.25 µl
for each analysis run	reaction vial	Master-Mix	water deionized
negative control NTC	colourless	12.75 µl	12.25 µl
positive control PTC	red	12.75 µl	12.25 µl

Tab. 1: PCRFast[®] Norovirus pipetting steps

–close the reaction vials (centrifuging recommended) and place them into the PCR Thermocycler

7.4. Instrument settings

Reporter: Fam and HEX
Quencher: non fluorescence
Reference dye: ROX (QuantiTect[™] Probe RT-PCR)

Note:

The QuantiTect[™] Probe RT-PCR master mix from Qiagen contains a high ROX level. Using realtime thermocycler which are able to excite ROX directly this results in extremely low dRn –values.

Perform amplification according to the following temperature / cycler profile:

Cycler profile

30 min	50 °C (122 °F)	reverse transcription
15 min	95 °C (203 °F)	activation of Taq polymerase
15 sec	95 °C (203 °F)	
60 sec	60 °C (140 °F)	45 cycles

Note:

Validation was performed using QuantiTect™ Probe RT-PCR MasterMix from Qiagen and the thermocycler specified in section 3. The specified cycler profile may need to be adjusted to the respective instrument and the MasterMix.

7.5. Detection

- assessment of the amplification in realtime is done via the amplification plot
- in principle the evaluation should be done using the software of the realtime thermocycler. Are the results not satisfactory using the instrument presettings, the “threshold” and the “baseline” settings could be change

8. Evaluation

The following table summarizes the evaluation. With a positive sample the final point of the fluorescence curve lies clearly above the threshold

8.1. Evaluation matrix

sample	reverse transcription control RTC	result
no amplification	amplification	sample negative
amplification	amplification	sample positive
no amplification	no amplification	inhibition *

negative control NTC	positive control PTC	result
no amplification		no contamination of master mix
	amplification	master mix is functional

* dilute extracted sample RNA once more and amplify again

Tab. 2: Evaluation of PCRFast® Norovirus

8.2.2. Examples of amplification plot evaluation


Fig. 1: sample negative; sample without amplification plot


Fig. 2: sample positive


Fig. 3: reaction inhibition, further diluting the isolated DNA

9. Sensitivity

The limit of detection is < 20 copies per reaction. Regarding norovirus in food the limit of detection is depending on the extraction method and efficiency of the reverse transcription.

10. Specificity

PCRFast[®] Norovirus is specific for norovirus genotype I and genotype II. The following species have been tested for cross - reactivity with at least 2500 copies each:

species		species		species	
Norovirus GG1	+	Norovirus GG2	+	Escherichia coli	-
Yersinia enterocolitica	-	Legionella pneumophilia	-	Campylobacter jejuni	-
Staphylococcus aureus	-	Bacillus cereus	-	Salmonellen	-
Clostridium perfringens	-	VTEC stx1	-	VTEC stx2	-
Legionella erythra	-	Shigella flexneri	-	Listeria monocytogenes	-
Clostridium difficile	-	Legionella pneumophilia	-	Helicobacter pylori	-

Tab. 3: Specificity of PCRFast[®] Norovirus

+ amplification

- no amplification

Should you have any questions on how to conduct this test, on sample preparations or on PCR analysis in general, please contact the Competence Centre for DNA Analysis of ifp, Institut für Produktqualität.

info@produktqualitaet.com

Tel. 0049 (0) 30 76 68 60 - 0

Fax 0049 (0) 30 76 68 60 - 50

You will also find more information on conducting PCRFast[®] at www.produktqualitaet.com.

These specifications are based on our current state of knowledge and are designed to provide information about our products and their potential applications. They are not intended to guarantee certain product properties or suitability for a specific application. ifp does not accept any liability except for the standard quality of the reagents. Replacement will be provided for any defective products. ifp is not liable for any other claims concerning direct or indirect damage or costs arising from the use of the products.

11. Literatur

Official collection of analysis methods after §64 LFGB

L 00.00 – 112 „Method for qualitative detection of Norovirus onto firm and smooth surfaces from food by Real-time PCR“

Dreier J., Störmer M., Mäde D., Burkhardt S., Kleesiek K. “Enhanced Reverse Transcription-PCR Assay for Detection of Norovirus Genogroup I” J Clin Microbiol (2006), **44 (8)**, 2714 - 2720.

Hoehne M., Schreier E. “Detection of Norovirus genogroup I and II by multiplex real-time RT-PCR using a 3'-minor groove binder-DNA probe” BMC Infectious Diseases (2006), **6 (69)**.
Microbiol., 73 (15), 2007, S. 4769-4775

ISO 22174

Microbiology of food and animal feeding stuffs - Polymerase chain reaction (PCR) for the detection of food - borne pathogens - General requirements and definitions

ISO 20837

Microbiology of food and animal feeding stuffs - Polymerase chain reaction (PCR) for the detection of food - borne pathogens - Requirements for sample preparation and qualitative detection

ISO 20838

Microbiology of food and animal feeding stuffs - Polymerase chain reaction (PCR) for the detection of food - borne pathogens - Requirements for amplification and detection for qualitative methods

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is a registered trademark of the ifp Institut für Produktqualität GmbH. ifp also offers contract analysis.