PCRFast[®] Clostridium perfringens

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22/07/2008

PCRFast[®] Clostridium perfringens

Realtime (SYBR[®] Green) and gel detection

PCR - test for the detection of Clostridium perfringens

Type No.: IF / MG1008

In-vitro test

Storage: 2 - 8 °C (35.6 - 46.4 °F)

Brief information

Simple, molecular biological test (PCR) for detecting Clostridium perfringens after enrichment from foodstuffs and animal feed in realtime (SYBR[®]Green) and / or gel detection. The kit detects the phospholipase C gene with an amplicon length of 283 bp. It could be also used for the identification of suspected colonies. With this kit 96 reactions could be done. Every reaction vial contains a specific primer pair. 48 reaction vials (red marking) contain additionally Clostridium perfringens DNA for investigating possible inhibiting effects (ITC) and for melting curve analysis.

Usage

For the general usage of the kit see "PCRFast[®] - Clostridium perfringens". Please read the instruction carefully.

Cycler profile

10 min	95 °C (203 °F)	45 cycles
15 sec	95 °C (203 °F)	
60 sec	56 °C (132.8 °F)	
30 sec	72 °C (161.6 °F)	

Sensitivity

The limit of detection is < 10 copies.

Specificity

PCRFast[®] Clostridium perfringens is 100 % specific for Clostridium perfringens. The following species have been tested with at least 2500 copies for cross reactivity:

Species		Species		Species	
Clostridium perfringens	+	Yersinia enterocolitica	-	Campylobacter jejuni	-
Listeria monocytogenes	-	Salmonellae	-	Salmonellae	-
Enterobacteriaceae STX1	-	Enterobacteriaceae STX2	-	Legionella erythra	-
Staphylococcus aureus	-	Bacillus cereus	-	Escherichia coli	-
Shigella flexneri	-				

Tab. 1: Specificity PCRFast[®] Clostridium perfringens

+ : band with 283bp

- : no band with 283bp

Official collection of analysis methods after § 64 LFGB

L 00.00 45 „General method specific requirements for the detection of microorganisms using the Polymerase Chain Reaction (PCR) in foodstuff “

Certificate of Quality

PCRFast[®] Clostridium perfringens,
complies specification.

Lot: **TCLO_11111**
quality assurance: A.Schuhr

22/07/2008