# PCRFast<sup>®</sup> Legionella pneumophila

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22/07/2008

# PCRFast<sup>®</sup> Legionella pneumophila

Realtime (SYBR<sup>®</sup> Green) and gel detection

## PCR - test for the detection of Legionella pneumophila

Type No.: IF / MG1021

In-vitro test

Storage: 2 - 8 °C (35.6 - 46.4 °F)

### Brief information

Simple, molecular biological test (PCR) for detecting Legionella pneumophila in realtime (SYBR<sup>®</sup> Green) and / or gel detection in water (16S rRNA gene with an amplicon size of 160 bp). With this kit 96 reactions could be done. Every reaction vial contains a specific primer pair. 48 reaction vials (red marking) contain additionally Legionella pneumophila DNA for investigating possible inhibiting effects (ITC) and for melting curve analysis.

### Usage

For the general usage of the kit see "PCRFast<sup>®</sup> - Legionella pneumophila ". Please read the instruction carefully.

### Cycler profile

10 min	95 °C (203 °F)	
15 sec	95 °C (203 °F)	
60 sec	72 °C (161.6 °F)	<b>45 cycles</b>

### Sensitivity

The limit of detection is 1 Legionella pneumophila in 1 ml water sample filtrating 100 ml.

### Specificity

PCRFast<sup>®</sup> Legionella pneumophila is 100 % specific for legionella. The following species have been tested with at least 2500 copies for cross reactivity:

Species		Species		Species	
Legionella pneumophila	+	Legionella erythra	-	Escherichia coli	-
Salmonella thyphimurium	-	Listeria monocytogenes	-	Campylobacter jejuni	-
Staphylococcus aureus	-	Bacillus cereus	-	Yersinia enterocolitica	-
Legionella pneumophila	-	Shigella flexneri	-		

Tab. 1: Specificity PCRFast<sup>®</sup> Legionella pneumophila

+ : band with 160 bp

- : no band with 160 bp

### ISO 22174

Microbiology of food and animal feeding stuffs - Polymerase chain reaction (PCR) for the detection of food - borne pathogens - General requirements and definitions

### Certificate of Quality

PCRFast<sup>®</sup> Legionella pneumophila ,  
complies specification.

Lot: TLEP\_11111  
quality assurance: A.Schuhr

22/07/2008