# PCRFast<sup>®</sup> Legionella pneumophila

Realtime (Sonde)

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“The Pneumoplex Assays, a Multiplex PCR-Enzyme Hybridization Assay That Allows Simultaneous Detection of Five Organisms, M. pneumoniae, C.pneumoniae, Legionella pneumophila, Legionella micdadei, and B. pertussis, and Its Real-Time Counterpart” J. Clin. Microbiol. Legionella (2005) 43, No. 2, p. 565–571

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# PCRFast<sup>®</sup> Legionella pneumophila

Realtime (probe)

## PCR - test for the detection of Legionella pneumophila

Type No.: IF / MR1021

In-vitro test

Storage: 2 - 8 °C (35.6 - 46.4 °F)

### Brief information

Easy-to-use molecular biological test (PCR) for the detection of Legionella pneumophila (probe with **Fam** reporter and non fluorescence quencher) in water (detection of the 16S rRNA gene) . The test kit has 96 reaction vials, each reaction vial contains a specific primer pair, the probe and an internal amplification control (probe with **HEX** reporter and non fluorescence quencher) for investigating possible inhibiting effects (ITC). 8 reaction vials (red marking) contain additionally Legionella pneumophila DNA for the PCR positive controls.

### Usage

For the general usage of the kit see the manual "PCRFast<sup>®</sup> - Legionella pneumophila". Please read the instruction carefully.

### Cycler profile

10 min	95 °C (203 °F)	Activation of polymerase
15 sec	95 °C (203 °F)	
30 sec	60 °C (140 °F)	<b>45 cycles</b>
30 sec	72 °C (161.6 °F)	

### Sensitivity

Limit of detection: 1 Legionella pneumophila in 100 ml sample volume.

### Specifity

PCRFast<sup>®</sup> Legionella pneumophila is 100 % specific for Legionella pneumophila. The following species have each been tested with at least 2500 copies for cross reactivity:

species		species		species	
Legionella pneumophila	+	Legionella erythra	-	Escherichia coli	-
Salmonella thyphimurium	-	Listeria monocytogenes	-	Campylobacter jejuni	-
Staphylococcus aureus	-	Bacillus cereus	-	Yersinia enterocolitica	-
Clostridium perfringens	-	Shigella flexneri	-		

Tab. 1: Specifity PCRFast<sup>®</sup> Legionella pneumophila

+ : amplification

- : no amplification

**Khanna M., Fan J., Pehler-Harrington K., Waters C., Douglass P., Stallock J., Kehl S., Henrickson K.J.**

"The Pneumoplex Assays, a Multiplex PCR-Enzyme Hybridization Assay That Allows Simultaneous Detection of Five Organisms, M. pneumoniae, C.pneumoniae, Legionella pneumophila, Legionella micdadei, and B. pertussis, and Its Real-Time Counterpart" J. Clin. Microbiol. Legionella (2005) 43, No. 2, p. 565-571

### ISO 22174

Microbiology of food and animal feeding stuffs - Polymerase chain reaction (PCR) for the detection of food - borne pathogens - General requirements

### Certificate of Quality

PCRFast<sup>®</sup> Legionella pneumophila,  
complies specification.

Lot: **RLEP\_11111**  
quality assurance: A.Schuhr

04/03/2009