



PURIFICATION OF FLUOROUS RADIOTRACERS

APPLICATION

AFFINIMIP® SPE ¹⁸F - Aromatic Nucleophilic Substitution is a selective SPE cartridge used to purify rapidly and perfectly a ¹⁸F radiotracer synthesized from an aromatic nucleophilic substitution. The cartridge retains the ammonium-based precursor as well as the phenol- and dialkylamino-based impurities.



RESULTS



Chromatograms obtained before (red) and after (E1 : blue, E2 : green) **AFFINIMIP® SPE ¹⁸F Aromatic Nucleophilic Substitution Clean-up**.

Recovery of more than 95 % of the fluororous radiotracer was obtained without any contamination of the other compounds.



PROTOCOL OF PURIFICATION

Sample preparation

A mixture of Ethyl 4-hydroxybenzoate (99 µg), Ethyl 4-dimethylaminobenzoate (928 µg), Ethyl 4-fluorobenzoate (10 µg) and Ethyl 4-trimethylammoniumbenzoate iodide (102 µg) in 80-20 Water-Acetonitrile (5 mL) is prepared as the loading solution.

Purification with a reversible **AFFINIMIP® SPE 18F Aromatic Nucleophilic Substitution** cartridge

Equilibration

- 5mL Acetonitrile
- 2mL Water

Loading

- 5mL of loading solution

Washing of interferences

- 5mL Water- Acetonitrile (80/20)

Drying 30 seconds

Elution (E)

- Elution of Ethyl 4-fluorobenzoate with 3mL of ACN (E1)
- Extra-elution of 1mL until dryness (E2)

HPLC Method with UV detection

Column: Hypersil Gold C18 column 50mm x 2.1mm

Mobile phase: gradient profile

| Time (min) | % water | % ACN |
|------------|---------|-------|
| 0 | 100 | 0 |
| 3 | 400 | 0 |
| 15 | 70 | 30 |
| 32 | 70 | 30 |
| 33 | 100 | 0 |
| 43 | 100 | 0 |

Flow rate: 0.2mL/min

UV detection: 235nm

Injection volume: 10µL.

Catalog number:

RP100-01 for 10 cartridges compatible with automates