## foodproof® STEC Screening LyoKit (-5' Nuclease-) [R 602 11-1] [R 602 11-2] Ver.1, Sep./2013

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**References**

1. Karch H, Tarr P, Bielaszewska M (2005). "Enterohaemorrhagic Escherichia coli in human medicine.". Int J Med Microbiol 295 (6-7): 405-18.
2. ISO/TS 13136:2012 "Microbiology of food and animal feed –Real-time polymerase chain reaction (PCR)–based method for the detection of food-borne pathogens – Horizontal method for the detection of Shiga toxin-producing Escherichia coli (STEC) and the determination of O157, O111, O26, O103 and O145 serogroups".
3. FSIS Microbiology Laboratory Guidebook (MLG) 5B, "Detection and Isolation of Non-O157 Shiga-Toxin Producing Escherichia coli (STEC) from Meat Products", 2011.

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**5. Supplementary Information**

**5.1 Ordering Information**

BIOTECON Diagnostics is offering a broad range of reagents and services. For a complete overview and for more information, please visit our website at [www.bc-diagnostics.com](http://www.bc-diagnostics.com).

**5.2 License**

**License Notice**

NOTICE TO PURCHASER: LIMITED LICENSE

Use of this product is covered by one or more of the following US patents and corresponding patent claims outside the US:5,804,375, 5,538,848, 5,723,591, 5,876,930, 6,030,787, 6,258,569. The purchase of this product includes a limited, non-transferable immunity from suit under the foregoing patent claims for using only this amount of product solely in Food Testing Applications and Genetically Modified Organism (GMO) Testing Applications, including reporting results of purchaser's activities for a fee or other commercial consideration, and also for the purchaser's own internal research. No right under any other patent claim is conveyed expressly, by implication, or by estoppel. Further information on purchasing licenses may be obtained from the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404, USA.

The purchase price of this product includes limited, nontransferable rights under U.S. Patent No. 7,687,247 owned by Life Technologies Corporation to use only this amount of the product to practice the claims in said patent solely for activities of the purchaser for bioburden testing, environmental testing, food testing, or testing for genetically modified organisms (GMO) in accordance with the instructions for use accompanying this product. No other rights are conveyed, including no right to use this product for *in vitro* diagnostic, therapeutic, or prophylactic purposes. Further information on purchasing licenses under the above patent may be obtained by contacting the Licensing Department, Life Technologies Corporation, 5791 Van Allen Way, Carlsbad, CA 92008. Email: [outlicensing@lifetech.com](mailto:outlicensing@lifetech.com).

**5.3 Trademarks**

foodproof® is a trademark of BIOTECON Diagnostics GmbH.

Other brand or product names are trademarks of their respective holders.

**5.4 Contact and Support**

If you have questions about this or any other product of BIOTECON Diagnostics, please contact our Technical Support staff (for details see [www.bc-diagnostics.com](http://www.bc-diagnostics.com)). Our scientists commit themselves to providing rapid and effective help. We also want you to contact us if you have suggestions for enhancing our product performance or using our products in new or specialized ways. Such customer information has repeatedly proven invaluable to us and the worldwide research community.

**6. Change Index**

Version 1, September 2013:

First version of the package insert.

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