

# DETERMINATION OF PATULIN IN TOMATO KETCHUP AND TOMATO POWDER

## PROTOCOL OF PURIFICATION

### Sample preparation

#### Preparation OF TOMATO KETCHUP

10g tomato ketchup and 10mL water are mixed with 150µL pectinase solution and left overnight at RT before a filtration with filter 0.2µm to obtain the loading solution.

#### Preparation OF TOMATO POWDER

10g tomato ketchup and 20mL water are mixed. 10g of the mixture, 10mL water and 150µL pectinase solution are left overnight at RT before a centrifugation at 4500rpm during 5 min. Then the mixture is filtered with filter 0.2µm to obtain the loading solution.

### Purification with a 3mL/100mg AFFINIMIP® SPE Patulin cartridge

#### Equilibration

- 2mL Acetonitrile
- 1mL Water

#### Loading

- 5mL of loading solution from tomato ketchup or 2mL from tomato powder

#### Washing of interferents (W1)

- 4mL Water-1% Acetic Acid
- 4mL Water

#### Drying by applying vacuum 10 seconds

#### Washing of interferents (W2)

- 500µL Diethyl Ether

#### Elution (E)

- 2mL Ethyl Acetate

The elution fraction was then evaporated and dissolved in water containing 0.1% acetic acid before HPLC analysis.

### HPLC Method


Column: Atlantis T3 column, 150mm x 2.1mm  
Mobile phase: gradient profile

Time (min)	% water	% ACN
0	98	2
20	98	2
21	50	50
25	50	50
26	98	2

Flow rate: 0.2mL/min  
Detection: UV - 276nm  
Injection volume: 100µL.


## RESULTS

### TOMATO KETCHUP



Chromatograms obtained after AFFINIMIP® SPE Patulin Clean-up of TOMATO KETCHUP spiked at 40µg/kg with Patulin (red) or not spiked (light blue).

### TOMATO POWDER



Chromatograms obtained after AFFINIMIP® SPE Patulin Clean-up of TOMATO POWDER spiked at 36µg/kg with Patulin (red) or not spiked (light blue).

### Catalog number:

#### 3mL-100mg sorbent

- FS102-02 for 25 cartridges
- FS102-03 for 50 cartridges
- FS102-02K for a kit of 25 cartridges + 50mL Pectinase
- FS102-03K for a kit of 50 cartridges + 50mL Pectinase
- REA-001-50mL for 50mL Pectinase solution