PCRFast[®] VTEC

**Geldetektion
gel detection**

IF/MG1004

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Kurzinformation

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11. Literatur

Herbert Schmidt, Caren Geitz, Phillip I. Tarr, Matthias Frosch, and Helge Karch, "Non-O157:H7 Pathogenic Shiga Toxin–Producing *Escherichia coli*: Phenotypic and Genetic Profiling of Virulence Traits and Evidence for Clonality", *The Journal of Infectious Diseases* 1999;179:115–23

Beutin, L., „Emerging Enterohaemorrhagic *Escherichia coli*, Causes und Effects of the Rise of a Human Pathogen“. *Journal Vet. Med. B* 53, Blackwell Verlag, Berlin, 2006, S. 299-303

Beutin, L.; Miko, A.; Krause, G.; Pries, K.; Haby, S.; Steege, K.; Albrecht, N. „Identification of human pathogenic strains of Shiga toxin-producing *Escherichia coli* from food by a combination of serotyping and molecular typing of Shiga toxin genes“. *Appl. Environ Microbiol.*, 73 (15), 2007, S. 4769-4775

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PCRFast[®] VTEC gel detection

Brief information

Simple, molecular biological test (polymerase chain reaction - PCR) for the qualitative detection of shigatoxin producing Escherichia coli (VTEC stx1 and stx2) in food with gel detection. With this kit 96 reactions could be done. Every reaction vial contains a specific primer pair.

48 reaction vials (red marking) contain additionally VTEC stx1 and VTEC stx2 DNA for investigating possible inhibition effects (ITC).

Time required:	enrichment.....	approx. 18 h
(10 samples)	DNA - Extraction.....	45 min
	PCR setup.....	15 min
	PCR.....	2.0 h
	gelelectrophoresis.....	15 - 30 min

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PCRFast[®]

is a registered trademark of ifp Institut für Produktqualität GmbH.
ifp also offers contract analyses.

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1. Principle of the test

PCRFast[®] VTEC is a simple, molecular biological test (PCR) for the qualitative detection of shigatoxin producing Escherichia coli (VTEC stx1 and stx2) in foodstuff. The analytical procedure described here complies with the international standards (e.g. ISO/DIS 20838) for PCR analysis.

All the reaction vials contain the optimum amount of the specific primer required for the PCR reaction.

The **reaction vials that are not marked (colourless)** are used for detecting **VTEC stx1 and VTEC stx2 DNA in the extracted sample** and for the **negative controls NTC** (checking the MasterMix for contamination).

The **marked reaction vials (red)** contain VTEC DNA as well as the primers. These are used for the **inhibition controls ITC** (checking the extracted sample DNA for inhibitors). The marked reaction vials are also used for the **positive controls PTC** (checking the functionality of MasterMix).

Place **12.5 µl of double concentrated MasterMix** into the reaction vial and add **12.5 µl of the DNA sample extract**. The MasterMix already contains polymerase, nucleotides and magnesium chloride in optimum concentrations. The target sequence is then amplified in a PCR thermocycler and the resulting amplicate (stx1: 282 bp / stx2: 584 bp) made visible in agarose gel electrophoresis using ethidium bromide.

2. Package contents

The test kit contains:

- 6 x **strips (colourless)** with 8 x 0.2 ml reaction vials each (48 reaction vials), coated with specific primers (for sample determinations and for the negative controls NTC)
- 6 x **strips (red)** with 8 x 0.2 ml reaction vials each (48 reaction vials), coated with specific primers and additionally VTEC stx1 and VTEC stx2 DNA (for the inhibition controls ITC and positive controls PTC)

3. Additionally required instruments and reagents

Instruments

- 37 °C (98.6 °F) incubator
- vortexer
- biocentrifuge, min. 14,000 x g (for 1.5 ml and 2.0 ml reaction vials)
- PCR thermocycler:**
 - gel, e.g. Eppendorf MasterCycler with heating block for 0.2 ml reaction vials
- gel chamber with power pack and gel documentation system
- microliter pipettes 2 - 20 µl, 20 - 200 µl, 100 - 1,000 µl, e.g. Gilson Pipetman P with filter tips
- 1.5 ml or 2.0 ml reaction vials

Reagents

Dilution buffer

- modified Tryptose Soy Bouillon with Novobiocin (mTSNB) according to LFGB § 64 L 00.00 - 92

DNA extraction

- ethanol >95 % (denatured)
- 0.1 x TE - buffer (TRIS 1 mmol / l, EDTA 0.1 mmol / l)
- water deionised, DNA free
- DNA-extraction kit, e.g. PCRFast[®] Microbe Extraction

PCR

for gel detection

double concentrated MasterMix (Hotstart recommended)

z. B. AmpliTaq Gold[®] PCR MasterMix, Applied Biosystems Nr. 4318739; Qiagen HotStarTaq MasterMix Kit Qiagen Nr. 203443; Brilliant II QPCR MasterMix, Stratagene Nr. 600804

Detection

- agarose
- ethidium bromide
- buffer for gel fabrication, gel overlaying, and electrophoresis running buffer: 1 x TAE - buffer (Tris 40 mmol / l, acetic acid 20 mmol / l, EDTA 1 mmol / l)
- gel loading buffer, e.g. 5 - fold (sucrose 5%, Orange G 0.25%)
- DNA length marker, range 50 - 1000 bp

4. Precautions

- all tasks should be performed taking all precautions common in labs, and all PCR tasks should be performed in accordance with the CEN / ISO recommendations
- to avoid carryovers, perform all work wearing protective gloves and use filter tips
- perform sample preparation, PCR setup and detection in separate rooms
- ethidium bromide is damaging to DNA and should therefore be handled in compliance with the regulations
- UV light is damaging to DNA, so take care when using the UV transilluminators

5. Storage instructions

Store the test kit and the reaction vials at 2 - 8 °C (35.6 - 46.4 °F).

6. Sample enrichment and DNA isolation

6.1. Enrichment

Foodstuff samples are enriched according to § 64 of the German Foodstuffs and Animal Feed Code (LFGB).

6.1.1. Enrichment

For enrichment of solid foodstuff samples are weighed in a sterile stomacher bag, diluted at a ratio of 1 : 10 (w / v) with dilution buffer (see 3.) and homogenised (e.g. 10 g sample + 90 ml dilution buffer). Incubated at 37 °C (98.6 °F) over night.

6.2. DNA isolation

For DNA extraction, take 1 ml from 6.1.1. and prepare the DNA accordingly to manual of the extraction kit manufacturer (e.g. the protocol of the PCRFast[®] Microbe Extraction kit).

6.3. Note

Pyrolyse is not advisable after enrichment in lauryl - sulfat broth.

6.4. Confirming suspect colonies

To confirm suspect colonies on solid culture medium, it suffices to slurry one loop of this colony in 0.2 ml of deionized water and heat it at 95 °C (203 °F) for 10 min. The sample is then cooled and centrifuged at 14,000 x g for 5 min. The supernatant can be used directly for PCRFast®.

7. PCR setup

7.1. Preparations

- remove strip from the plastic bag and separate required number of reaction vials
- place the remaining strips / reaction vials and the desiccant back into the bag, seal bag and store at 2 - 8 °C (35.6 - 46.4 °F)
- have the required amount of MasterMix at hand (12.5 µl per reaction)

7.2. Procedure

- pipette the following volumes into the reaction vials:

per sample	reaction vial	MasterMix	DNA from sample
	colourless	12.5 µl	12.5 µl
inhibition control ITC	red	12.5 µl	12.5 µl
per analysis run	reaction vial	Master-Mix	0.1 x TE buffer
negative control NTC	colourless	12,5 µl	12.5 µl
positive control PTC	red	12.5 µl	12.5 µl

Tab. 1: PCRFast® VTEC pipetting steps

- close the reaction vials (centrifuging recommended) and place them into the PCR thermocycler. Perform amplification according to the following temperature / cycler profile:

Cycler profile

15 min	95 °C (203 °F)	
15 sec	95 °C (203 °F)	
30 sec	60 °C (140 °F)	45 cycles
45 sec	72 °C (161.6 °F)	

Note:

Validation was performed using AmpliTaq Gold® PCR MasterMix by Applied Biosystems and the thermocycler specified in section 4. The specified cycler profile may need to be adjusted to the respective instrument and the MasterMix.

7.3. Detection

7.3.1. Gel detection

- we recommend electrophoresis in 2.0 - 2.5 % agarose and staining the amplificate with ethidium bromide
- perform gel electrophoresis according to the manufacturer's instructions

Example:

- fabricate agarose gel 2.0 % - 2.5 % (dissolve 1.0 - 1.25 g agarose and 2 - 4 µl ethidium bromide 10 mg / ml in 50 ml of 1 x TAE buffer by heating and pour gel
- overlay agarose gel with 1 x TAE buffer
- apply PCR amplificate with gel loading buffer to the agarose gel, e.g. add 6 µl of the 5 x loading buffer to the PCR amplificate (25 µl), mix and apply 10 µl
- apply length markers
- perform electrophoresis at 3 - 6 V / cm (depending on electrophoresis system) for 15 - 30 min
- visualize DNA on a transilluminator and document results

8. Evaluation

8.1. Gel detection

- The following table summarizes the evaluation. With a positive sample, the sample band of one toxin or the sample bands of both toxin and the bands of the positive control are at the same level.

8.1.1. Evaluation matrix for gel detection


Sample	Inhibition control ITC	negative control NTC	positive control PTC	Result
██████	██████			sample positive
no band	██████			sample negative
no band	no band			inhibition*
		no band		MasterMix not contaminateded
			██████	MasterMix functional

* dilute extracted sample DNA once more and amplify again; ██████ = band

Tab. 2: Evaluation of PCRFast® VTEC gel detection

8.1.2. Example of a typical agarose gel image

– The following figure shows an exemplary agarose gel image.

sample 1 to sample 2: VTEC stx2 positive
 sample 3: VTEC stx1 und stx2 positive (equal to positive control)
 sample 4 to sample 5: VTEC stx 1 positive
 NTC: negativ control
 M: marker (50 bp, 150 bp, 300 bp, 500pb, 766 bp)

Fig. 1: Example of a PCRFast® VTEC agarose gel image

9. Sensitivity

The Limit of detection is < 10 copies per reaction. Regarding shigatoxin producing Escherichia coli (VTEC stx 1 / 2) in food the detection limit is less than 1 shigatoxin producing Escherichi coli in 25 g sample.

10. Specificity

PCRFast[®] VTEC is specific for Shigatoxin producing Escherichia coli (VTEC stx1 and stx2). The following species have been tested for cross - reactivity with at least 2500 copies each:

Species		Species		Species	
VTEC stx1	+	VTEC stx2	+	Escherichia coli	-
Yersinia enterocolitica	-	Legionella pneumophila	-	Campylobacter jejuni	-
Staphylococcus aureus	-	Bacillus cereus	-	Salmonellae	-
Clostridium perfringens	-	Listeria monocytogenes	-	Legionella erythra	-

Tab. 4: Specificity of PCRFast[®] VTEC

- + band with 282 / 584 bp
- + no band with 282 / 584 bp

Should you have any questions on how to conduct this test, on sample preparations or on PCR analysis in general, please contact the Competence Centre for DNA Analysis of ifp, Institut für Produktqualität.

info@produktqualitaet.com

Tel. 0049 (0) 30 76 68 60 - 0

Fax 0049 (0) 30 76 68 60 - 50

You will also find more information on conducting PCRFast[®] at www.produktqualitaet.com.

These specifications are based on our current state of knowledge and are designed to provide information about our products and their potential applications. They are not intended to guarantee certain product properties or suitability for a specific application. ifp does not accept any liability except for the standard quality of the reagents. Replacement will be provided for any defective products. ifp is not liable for any other claims concerning direct or indirect damage or costs arising from the use of the products.

11. Literatur

Herbert Schmidt, Caren Geitz, Phillip I. Tarr, Matthias Frosch, and Helge Karch, “Non-O157:H7 Pathogenic Shiga Toxin–Producing *Escherichia coli*: Phenotypic and Genetic Profiling of Virulence Traits and Evidence for Clonality”, *The Journal of Infectious Diseases* 1999;179:115–23

Beutin, L., „Emerging Enterohaemorrhagic *Escherichia coli*, Causes und Effects of the Rise of a Human Pathogen“. *Journal Vet. Med. B* 53, Blackwell Verlag, Berlin, 2006, S. 299-303

Beutin, L.; Miko, A.; Krause, G.; Pries, K.; Haby, S.; Steege, K.; Albrecht, N. „Identification of human pathogenic strains of Shiga toxin-producing *Escherichia coli* from food by a combination of serotyping and molecular typing of Shiga toxin genes“. *Appl. Environ Microbiol.*, 73 (15), 2007, S. 4769-4775

ISO 22174

Microbiology of food and animal feeding stuffs - Polymerase chain reaction (PCR) for the detection of food - borne pathogens - General requirements and definitions

ISO 20837

Microbiology of food and animal feeding stuffs - Polymerase chain reaction (PCR) for the detection of food - borne pathogens - Requirements for sample preparation and qualitative detection

ISO 20838

Microbiology of food and animal feeding stuffs - Polymerase chain reaction (PCR) for the detection of food - borne pathogens - Requirements for amplification and detection for qualitative methods

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is a registered trademark of the ifp Institut für Produktqualität GmbH.