PCRFast[®]
Enterobacter sakazakii

Realtime (SYBR[®] Green) und / oder Gel
Realtime (SYBR[®] Green) and / or Gel

IF/MG1019

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Kurzinformation

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11. Literatur

Seo, K. H. et al.

“Rapid, Specific Detection of Enterobacter sakazakii in Infant Formula Using a Real – Time PCR Assay” Journal of Food Protection, 68 No.1, 59-63

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PCRFast[®] Enterobacter sakazakii Realtime (SYBR[®]Green) and / or gel detection

Brief information

Simple, molecular biological test (PCR) for detecting Enterobacter sakazakii in foodstuffs and animal feed in realtime (SYBR[®]Green) and / or gel detection (detection of the intergenic sequenz from rpsU - dnaG gene with 78 bp). The test kit has 96 reaction vials, each reaction vial contains a specific primer pair.

48 reaction vials (red marking) contain additionally Enterobacter sakazakii DNA for investigating possible inhibiting effects (ITC) and for melting curve analysis.

Time required:	enrichment.....	approx. 24 h
	extraction.....	25 min
(10 samples)	PCR setup.....	15 min
	PCR.....	2.0 h
	gelelectrophoresis.....	15 - 30 min

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PCRFast[®]
is a registered trademark of ifp Institut für Produktqualität GmbH.
ifp also offers contract analyses.

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1. Principle of the test

PCRFast[®] Enterobacter sakazakii is a simple, molecular biological test (PCR) for detecting Enterobacter sakazakii in foodstuffs and animal feed (detection of the intergenic sequence from rpsU - dnaG gene with 78 bp). The analytical procedure described here complies with the international standards (ISO) for PCR analysis and the § 64 of the German Foodstuffs and Animal Feed Code (LFGB).

All the reaction vials contain the optimum amount of the specific primer required for the PCR reaction.

The **reaction vials that are not marked (colourless)** are used for detecting **Enterobacter sakazakii DNA in the extracted sample** and for the **negative controls NTC** (checking the MasterMix for contamination).

The **marked reaction vials (red)** contain Enterobacter sakazakii DNA as well as the primers. These are used for the **inhibition controls ITC** (checking the extracted sample DNA for inhibitors) and for the melting curve analysis. The marked reaction vials are also used for the **positive controls PTC** (checking the functionality of MasterMix).

Place **12.5 µl of double concentrated MasterMix** into the reaction vial and add **12.5 µl of the DNA sample extract**. The MasterMix already contains polymerase, nucleotides and magnesium chloride in optimum concentrations. In realtime the MasterMix contains additionally SYBR[®]Green. The target sequence is then amplified in a PCR thermocycler and the resulting amplificate made visible in realtime (SYBR[®]Green) and / or agarose gel electrophoresis using ethidium bromide. With realtime (SYBR[®]Green), positive reactions need to be verified via melting curve analysis. To do so, the melting curve peaks of the sample are compared to those of the **inhibition control ITC**. The amplificate can additionally be verified by means of agarose gel electrophoresis via the base length.

2. Package contents

The test kit contains:

- 6 x **strips (colourless)** with eight 0.2 ml reaction vials each (48 reaction vials), coated with specific primers (for sample determinations and for the negative controls NTC)
- 6 x **strips (red)** with eight 0.2 ml reaction vials each (48 reaction vials), coated with specific primers and additionally *Enterobacter sakazakii* DNA (for the inhibition controls ITC and positive controls PTC)

3. Additionally required instruments and reagents

Lab material and instruments

- homogenizer (Stomacher)
- stomacher bag, sterile
- measuring cylinder, sterile
- erlenmeyer flask, sterile
- culture tube, sterile
- 37 °C (98.6 °F) incubator
- vortexer
- biocentrifuge, min. 14,000 x g (for 1.5 ml and 2.0 ml reaction vials)
- PCR thermocycler:**
 - gel, e.g. Eppendorf MasterCycler with heating block for 0.2 ml reaction vials
 - Realtime, e.g. Applied Biosystems 7500, STRATAGENE Mx 3005 P, Mx 3000 P with software for the melting curve analysis each
- gel chamber with power pack and gel documentation system
- microliter pipettes 2 - 20 µl, 20 - 200 µl, 100 - 1,000 µl, e.g. Gilson Pipetman P, with filter tips
- 1.5 ml or 2.0 ml reaction vials

Reagents

Pre - Enrichment

- sterile deionized water (prewarmed to 45 °C)

Enrichment

- Enterobacteriaceae enrichment broth (EE - broth; Becton Dickinson)

DNA extraction

- ethanol >95 % (denatured)
- 0.1 x TE - buffer (TRIS 1 mmol / l, EDTA 0.1 mmol / l)
- water deionised, DNA free
- DNA-extraction kit, e.g. PCRFast[®] Microbe Extraction

PCR

only for gel detection

double concentrated MasterMix (Hotstart recommended)

e.g. AmpliTaq Gold[®] PCR MasterMix, Applied Biosystems No. 4318739; Qiagen HotStarTaq MasterMix Kit Qiagen No. 203443; Brilliant II QPCR MasterMix, Stratagene No. 600804

for gel detection and / or Realtime (SYBR[®]Green)

double concentrated MasterMix (Hotstart recommended)

e.g. Power SYBR[®] Green PCR MasterMix, Applied Biosystems No. 4367659; QuantiTect SYBR[®]Green PCR Kit Qiagen No. 204143; Brilliant II SYBR[®] Green MasterMix Stratagene No. 600548

Detection

- agarose
- ethidium bromide
- buffer for gel fabrication, gel overlaying, and electrophoresis running buffer:
 - 1 x TAE - buffer (Tris 40 mmol / l, acetic acid 20 mmol / l, EDTA 1 mmol / l)
- gel loading buffer, e.g. 5 - fold (sucrose 5 %, Orange G 0.25 %)
- DNA length marker, range 50 - 1000 bp

4. Precautions

- all tasks should be performed taking all precautions common in labs, and all PCR tasks should be performed in accordance with the CEN / ISO recommendations
- to avoid carryovers, perform all work wearing protective gloves and use filter tips
- perform sample preparation, PCR setup and detection in separate rooms
- ethidium bromide is damaging to DNA and should therefore be handled in compliance with the regulations
- UV light is damaging to DNA, so take care when using the UV transilluminators

5. Storage instructions

Store the test kit and the reaction vials at 2 - 8 °C (35.6 - 46.4 °F).

6. Sample enrichment and DNA isolation

6.1. Pre - Enrichment

25 g (ml) of a foodstuff sample is weighed in a sterile stomacher bag, diluted at a ratio of 1 : 10 (w / v) with sterile deionized water (prewarmed) and incubated at 37 °C (98.6 °F) over night (e.g. 25 g sample + 225 ml water).

With each analyzing series one sample without matrix should be done as an extraction control. A positive extraction control indicates a contamination in the reagents.

6.1.1. Selective enrichment

10 ml of the pre - enrichment from 6.1. is transferred to 90ml Enterobacteriaceae – Enrichment – Broth and incubated at 37 °C (98.6 °F) 24 h.

6.2. DNA isolation

Thermal Extraction procedure

- 1) for DNA - extraction take 1 ml from the enrichment culture (6.1.2.) and centrifuge at 14,000 x g for 10 min, discard supernatant
- 2) resuspend the pellet in 200 µl 0,1 x TE – buffer, centrifuge at 14,000 x g for 10 min, discard supernatant
- 3) resuspend the pellet again in 200 µl 0,1 x TE – buffer
- 4) thermal lysis: heat the pellet for 15 min. at 95 °C (203 °F)
- 5) cool down the solution (e.g. 1-2 min at - 20 °C (- 4 °F))
- 6) mix the solution and centrifuge at 14,000 x g for 5 min
- 7) transfer the supernatant (contains the DNA) in a new reaction vial

After the extraction the supernatant has to be diluted with a ratio of 1 : 10 with 0.1 x TE – buffer and can be used for the PCRFast[®] reaction (e.g. 10 µl supernatant + 90 µl 0,1 x TE – buffer).

Moreover every conventional extraktionsystem can be used for the DNA – extraction (e.g. PCRFast[®] Microbe Extraction).

6.3. Confirming suspect colonies

To confirm suspect colonies on solid culture medium, it suffices to slurry one loop of this colony in 0.2 ml of deionized water and heat it at 95 °C (203 °F) for 10 min. The sample is then cooled and centrifuged at 14,000 x g for 5 min. The supernatant can be used directly for PCRFast[®].

7. PCR setup

7.1. Preparations

- remove strip from the plastic bag and separate required number of reaction vials
- place the remaining strips / reaction vials and the desiccant back into the bag, seal bag and store at 2 - 8 °C (35.6 - 46.4 °F)
- have the required amount of MasterMix at hand (12.5 µl per reaction)

7.2. Procedure

- pipette the following volumes into the reaction vials:

per sample	reaction vial	MasterMix	DNA from sample
	colourless	12.5 µl	12.5 µl
inhibition control ITC	red	12.5 µl	12.5 µl
per analysis run	reaction vial	Master-Mix	0.1 x TE buffer
negative control NTC	colourless	12,5 µl	12.5 µl
positive control PTC	red	12.5 µl	12.5 µl

Tab. 1: PCRFast[®] Enterobacter sakazakii pipetting steps

- close the reaction vials (centrifuging recommended) and place them into the PCR thermocycler. Perform amplification according to the following temperature / cycler profile:

Cycler profile

10 min	95 °C (203 °F)	
15 sec	95 °C (203 °F)	
60 sec	62 °C (143.6 °F)	45 cycles

Note:

Validation was performed using AmpliTaq Gold[®] PCR MasterMix (for the gel detection) and Power SYBR[®] Green PCR MasterMix (for realtime and / or gel detection) by Applied Biosystems and the thermocyclers specified in section 3. The specified cycler profile may need to be adjusted to the respective instrument and the MasterMix.

7.3. Detection

7.3.1. Gel detection

- we recommend electrophoresis in 2.0 - 2.5 % agarose and staining the amplificate with ethidium bromide
- perform gel electrophoresis according to the manufacturer's instructions

Example:

- fabricate agarose gel 2.0 % - 2.5 % (dissolve 1.0 - 1.25 g agarose and 2 - 4 µl ethidium bromide 10 mg / ml in 50 ml of 1 x TAE buffer by heating and pour gel
- overlay agarose gel with 1 x TAE buffer
- apply PCR amplificate with gel loading buffer to the agarose gel, e.g. add 6 µl of the 5 x loading buffer to the PCR amplificate (25 µl), mix and apply 10 µl
- apply length markers
- perform electrophoresis at 3 - 6 V / cm (depending on electrophoresis system) for 15 - 30 min
- visualize DNA on a transilluminator and document results

7.3.2. Realtime (SYBR[®] Green)

- assessment of the amplification in realtime (SYBR[®] Green) is done via the amplification plot and the melting curve analysis
- negative samples can be sorted out

- if the amplification plot shows an increase in fluorescence, the sample must be verified via the melting curve analysis
- to do so, compare the sample melting curve with the melting curve resulting from the inhibition control ITC
- with positive samples, the sample's melting curve peak must agree with the melting curve peak of the inhibition control ITC (temperature deviance not more than ± 0.5 °C (0.9 °F))
- in order to additionally verify a positive sample via the base length, we recommend performing an electrophoresis in 2.0 % - 2.5 % agarose gel and staining the amplification product with ethidium bromide

8. Evaluation

8.1. Gel detection

The following table summarizes the evaluation. With a positive sample, the sample band and the band of the positive control are at the same level.

8.1.1. Evaluation matrix for gel detection

Sample	Inhibition control ITC	negative control NTC	positive control PTC	Result
██████	██████			sample positive
no band	██████			sample negative
no band	no band			inhibition*
		no band		MasterMix not contaminateded
			██████	MasterMix functional


* dilute extracted sample DNA once more and amplify again

██████ = band

Tab. 2: Evaluation of PCRFast[®] Enterobacter sakazakii gel detection

8.1.2. Example of a typical agarose gel image

– The following figure shows an exemplary agarose gel image.

upper lanes 1 and 2
(reaction vials colourless, containing specific primers)
sample 1: positive
sample 2: negative

lower lanes 3 and 4
(reaction vials red, containing specific primers and specific *Enterobacter sakazakii* DNA)
sample 1: no inhibition
sample 2: no inhibition

Fig. 1: Example of a PCRFast[®] *E.sakazakii* agarose gel image

8.2. Realtime (SYBR[®] Green)

8.2.1. Realtime (SYBR[®] Green) evaluation matrix

Sample	Inhibition control (ITC)	Result
no amplification	amplification + melting curve	sample negative
amplification + melting curve, but different melting temperature than inhibition control ITC (≥ 0.5 °C)	amplification + melting curve	sample negative
amplification + melting curve with same melting temperature as inhibition control ITC (± 0.5 °C)	amplification + melting curve	sample positive

Negative control NTC	Positive control PTC	Result
no amplification	no amplification	inhibition*
no amplification		MasterMix not contaminated
	amplification + melting curve	MasterMix functional

* dilute extracted sample DNA once more and amplify again

Tab. 3: Evaluation of PCRFast® E.sakazakii (Realtime SYBR® Green)

8.2.2. Examples of amplification plot evaluation


Fig. 2: sample negative; sample without amplification plot


Fig. 3: sample positive, but result verification via melting curve analysis


Fig. 4: reaction inhibition

8.2.3. Examples of melting curve evaluation


Fig. 5: sample negative; sample without melting curve


Fig. 6: sample positive; sample and inhibition control ITC with same melting curve


Fig. 7: sample negative; sample and inhibition control ITC with different melting curves

9. Sensitivity

Limit of detection: less than 1 *Enterobacter sakazakii* in 25 g sample.

10. Specificity

PCRFast[®] *Enterobacter sakazakii* is specific for *Enterobacter sakazakii* and *Shigella*. The following species have each been tested for cross-reactivity with at least 2500 copies:

species		species		species	
<i>Escherichia coli</i>	-	<i>Enterobacter sakazakii</i>	+	<i>Legionella pneumophila</i>	
Salmonellen	-	<i>Listeria monocytogenes</i>	-	<i>Campylobacter jejuni</i>	-
<i>Staphylococcus aureus</i>	-	<i>Bacillus cereus</i>	-	<i>Yersinia enterocolitica</i>	-
<i>Clostridium perfringens</i>	-	<i>Legionella erythra</i>	-	<i>Shigella flexneri</i>	-
VTEC stx1	-	VTEC stx2	-		

Tab. 4: Specificity of PCRFast[®] *Enterobacter sakazakii*

+ band with 78 bp

- no band with 78 bp

Should you have any questions on how to conduct this test, on sample preparations or on PCR analysis in general, please contact the Competence Centre for DNA Analysis of ifp, Institut für Produktqualität.

info@produktqualitaet.com

Tel. 0049 (0) 30 76 68 60 - 0

Fax 0049 (0) 30 76 68 60 - 50

You will also find more information on conducting PCRFast[®] at www.produktqualitaet.com.

These specifications are based on our current state of knowledge and are designed to provide information about our products and their potential applications. They are not intended to guarantee certain product properties or suitability for a specific application. ifp does not accept any liability except for the standard quality of the reagents. Replacement will be provided for any defective products. ifp is not liable for any other claims concerning direct or indirect damage or costs arising from the use of the products.

11. References

Seo, K. H. et al.

“Rapid, Specific Detection of *Enterobacter sakazakii* in Infant Formula Using a Real – Time PCR Assay” *Journal of Food Protection*, 68 No.1, 59-63

ISO 22174

Microbiology of food and animal feeding stuffs - Polymerase chain reaction (PCR) for the detection of food - borne pathogens - General requirements and definitions

ISO 20837

Microbiology of food and animal feeding stuffs - Polymerase chain reaction (PCR) for the detection of food - borne pathogens - Requirements for sample preparation and qualitative detection

ISO 20838

Microbiology of food and animal feeding stuffs - Polymerase chain reaction (PCR) for the detection of food - borne pathogens - Requirements for amplification and detection for qualitative methods

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is a registered trademark of the ifp Institut für Produktqualität GmbH. ifp also offers contract analysis.