PCRFast[®]
Legionella pneumophila

Realtime (Sonde)
Realtime (probe)

IF/MR1021

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Kurzinformation

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“The Pneumoplex Assays, a Multiplex PCR-Enzyme Hybridization Assay That Allows Simultaneous Detection of Five Organisms, *M. pneumoniae*, *C.pneumoniae*, *Legionella pneumophila*, *Legionella micdadei*, and *B. pertussis*, and Its Real-Time Counterpart” *J. Clin. Microbiol.* (2005) 43, No. 2, p. 565–571

PCRFast[®] Legionella pneumophila Realtime (probe)

Brief information

Simple, molecular biological test (PCR) for detecting Legionella pneumophila in realtime (probe with Fam reporter and non fluorescence quencher) in water (detection of the 16S rRNA gene). The test kit has 96 reaction vials, each reaction vial contains a specific primer pair, the probe and an internal amplification control (probe with HEX reporter and non fluorescence quencher) for investigating possible inhibiting effects (ITC).

8 reaction vials (red marking) contain additionally Legionella pneumophila DNA for the PCR positive controls (PTC).

Time required:	enrichment.....	approx. 24 h
	extraction.....	25 min
(10 samples)	PCR setup.....	15 min
	PCR.....	1,5 h

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PCRFast[®]
is a registered trademark of ifp Institut für Produktqualität GmbH.
ifp also offers contract analyses.

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1. Principle of the test

PCRFast[®] Legionella pneumophila is a simple, molecular biological test (PCR) for detecting Legionella pneumophila in water (detection of the 16SrRNA gene). The analytical procedure described here complies with the international standards (ISO) for PCR analysis and the § 64 of the German Foodstuffs and Animal Feed Code (LFGB).

All reaction vials contain the for the PCR reaction necessary specific primers and probe in an optimum amount and an internal amplification control for investigating inhibitory effects.

The **88 colourless reaction vials** contain specific primers and probe (Fam reporter and non fluorescence quencher) and are used for the detection of **Legionella pneumophila DNA in the sample extract**, for the **negative control NTC** (check on MasterMix for contamination) and for the **extraction control ETC** (check on extraction for contamination).

The **8 red - coloured reaction vials** contain beside the primers, probe and ITC additionally Legionella pneumophila DNA. These reaction vials are used for the **positive control PTC** (check on functionality of MasterMix).

12.5 µl of double concentrated MasterMix and, subsequently, **12.5 µl of the extracted DNA sample** is added into the reaction vial. The MasterMix contains an optimum concentration of polymerase, nucleotides and magnesium chloride. The target sequence is then amplified in a PCR thermocycler. Is the target present the probe hybridises on the built up DNA fragments. Thereby the reporter (Fam/HEX) and the quencher (non fluorescence) of the probe will be seperated and the reporter radiates a fluorescence signal (Fam: 520 nm, HEX: 553nm). The intensity of the signal increases with the amount of synthesised DNA which could be measured and displayed. If the target sequence is not present no increase in the fluorescence is detectable.

2. Package contents

The test kit contains:

- 11 x **strips (colourless)** with 8 x 0.2 ml reaction vials each (88 reaction vials), coated with specific primers (for sample determinations, the negative controls NTC and for the extraction control ETC) and an internal amplification control.
- 1 x **strips (red)** with 8 x 0.2 ml reaction vials each (8 reaction vials), coated with specific primers, probe and additionally Legionella pneumophila DNA (for the positive controls PTC) and an internal amplification control.

3. Additionally required instruments and reagents

Lab material and instruments

- filtration equipment
- measuring cylinder, sterile
- Petri dishes Ø 60mm, sterile
- 37 °C (98.6 °F) and 70°C (158 °F) incubator
- vortexer
- biocentrifuge, min. 14,000 x g (for 1.5 ml and 2.0 ml reaction vials)
- PCR thermocycler:**
 - e.g. Applied Biosystems 7500, STRATAGENE Mx 3005 P, Mx 3000 P
- microliter pipettes 2 - 20 µl, 20 - 200 µl, 100 - 1,000 µl, e.g. Gilson Pipetman P, each with filter tips
- 1.5 ml or 2.0 ml reaction vials

Reagents

Filtration

- 0.45 µm nitrocellulose filter Ø 50mm (or appropriate material),
e.g. Sartorius Celulose Nitrate Filter Order no.: 11306 50 N

DNA extraction

- ethanol >95 % (denatured)
- 0.1 x TE - buffer (TRIS 1 mmol / l, EDTA 0.1 mmol / l)
- water deionised, DNA free
- DNA-extraction kit, e.g. PCRFast[®] Microbe Extraction

PCR

Realtime

double concentrated MasterMix (Hotstart recommended),
e.g. AmpliTaq Gold[®] PCR MasterMix, Applied Biosystems No. 4318739;
Brilliant[®] II QPCR MasterMix, Stratagene No. 600804; QuantiTect Probe
PCR Kit, Qiagen No. 204343

4. Precautions

- all tasks should be performed taking all precautions common in labs, and all PCR tasks should be performed in accordance with the CEN / ISO recommendations
- to avoid carryovers, perform all work wearing protective gloves and use filter tips
- perform sample preparation, PCR setup and detection in separate rooms

5. Storage instructions

Store the test kit and the reaction vials at 2 - 8 °C (35.6 - 46.4 °F).

6. Sample enrichment and DNA isolation

6.1. Filtration

100 ml of the water sample will be filtrated through a 0.45 µm nitrocellulose filter Ø 50 mm with a filtration equipment.

With each analyzing series one sample without matrix should be done as an extraction control. A positive extraction control indicates a contamination in the reagents.

6.2. DNA isolation

The recommendation for DNA – Extraction is depended on the PCRFast[®] Microbe extraction Kit. For the DNA-extraction place the filter with 1225 µl of lysis puffer (e.g. PCRFast[®] Microbe Extraction 600 µl M1 and M2 each and 25 µl Proteinase K) in one Petridish Ø 60mm. Incubate the dish for 15 minutes at 70 °C (158°F) with shaking. Transfer the lysate completely (approx. 600 µl – 800 µl) in a reaction tube and mix it with Ethanol in a ratio of 2 : 1 (e.g. 800 µl lysat + 400 µl Ethanol). Apply entire solution to the column, centrifuge at 11,000 x g for 1 min; discard flow-through. If it is not possible to load the

complete solution to the column repeat this step. Continue the extraction procedure according to the manual of the extraction kit manufacturer (e.g. PCRFast[®] Microbe Extraction). After the extraction the eluate has to be diluted with a ratio of 1 : 10 with 0.1 x TE buffer.

6.3. Confirming suspect colonies

To confirm suspect colonies on solid culture medium, it suffices to slurry one loop of this colony in 0.2 ml of deionized water and heat it at 95 °C (203 °F) for 10 min. The sample is then cooled and centrifuged at 14,000 x g for 5 min. The supernatant can be used directly for PCRFast[®].

7. PCR setup

7.1. Preparations

- remove strip from the plastic bag and separate required number of reaction vials
- place the remaining strips / reaction vials and the desiccant back into the bag, seal bag and store at 2 - 8 °C (35.6 - 46.4 °F)
- have the required amount of MasterMix at hand (12.5 µl per reaction)

7.2. Procedure

–pipette the following volumes into the reaction vials:

per sample	reaction vial	MasterMix	DNA from sample
sample	colourless	12.5 µl	12.5 µl
per analysis run	reaction vial	Master-Mix	0.1 x TE buffer
negative control NTC	colourless	12,5 µl	12.5 µl
positive control PTC	red	12.5 µl	12.5 µl

Tab. 1: PCRFast® Legionella pneumophila pipetting steps

–close the reaction vials (centrifuging recommended) and place them into the PCR thermocycler.

7.4. Instrument settings

Reporter: Fam and HEX
Quencher: non fluorescence
Reference dye: ROX (depending on used master mix)

Perform amplification according to the following temperature / cycler profile:

Cycler profile

10 min	95 °C (203 °F)	
15 sec	95 °C (203 °F)	
30 sec	60 °C (140 °F)	45 cycles
30 sec	72 °C (161.6 °F)	

Note:

Validation was performed using Brilliant® II QPCR Master Mix by Stratagene and the thermocyclers specified in section 3. The specified cycler profile may need to be adjusted to the respective instrument and the MasterMix.

7.5. Detection

- assessment of the amplification in realtime is done via the amplification plot
- in principle the evaluation should be done using the software of the realtime thermocycler. Are the results not satisfactory using the instrument presettings, the threshold and the baseline settings could be changed.

8. Evaluation

The following table summarizes the evaluation. With a positive sample the final point of the fluorescence curve lies clearly above the threshold

8.1. Evaluation matrix

Sample	Inhibition control (ITC)	Negative control NTC	Positive control PTC	Result
+	+			sample positive
no amplification	+			sample negative
no amplification	no amplification			inhibition*
		no amplification		MasterMix not contaminated
			+	MasterMix functional

* dilute extracted sample DNA once more and amplify again

Tab. 3: Evaluation of PCRFast® Legionella pneumophila

8.2.2. Examples of amplification plot evaluation


Fig. 1: sample negative; sample without amplification plot


Fig. 2: sample positive


Fig. 3: reaction inhibition, further diluting the isolated DNA

9. Sensitivity

Limit of detection: 1 *Legionella pneumophila* in 1ml water sample
filtrating 100 ml.

10. Specificity

PCRFast[®] *Legionella pneumophila* is 100% specific for *Legionella pneumophila*. The following species have each been tested for cross-reactivity with at least 2500 copies:

species		species		species	
<i>Legionella pneumophila</i>	+	<i>Legionella erythra</i>	-	<i>Escherichia coli</i>	-
<i>Salmonella thyphimurium</i>	-	<i>Listeria monocytogenes</i>	-	<i>Campylobacter jejuni</i>	-
<i>Staphylococcus aureus</i>	-	<i>Bacillus cereus</i>	-	<i>Yersinia enterocolitica</i>	-
<i>Clostridium perfringens</i>	-	<i>Shigella flexneri</i>	-		

Tab. 4: Specificity of PCRFast[®] *Legionella pneumophila*

+ amplification

+ no amplification

Should you have any questions on how to conduct this test, on sample preparations or on PCR analysis in general, please contact the Competence Centre for DNA Analysis of ifp, Institut für Produktqualität.

info@produktqualitaet.com

Tel. 0049 (0) 30 76 68 60 - 0

Fax 0049 (0) 30 76 68 60 - 50

You will also find more information on conducting PCRFast[®] at www.produktqualitaet.com.

These specifications are based on our current state of knowledge and are designed to provide information about our products and their potential applications. They are not intended to guarantee certain product properties or suitability for a specific application. ifp does not accept any liability except for the standard quality of the reagents. Replacement will be provided for any defective products. ifp is not liable for any other claims concerning direct or indirect damage or costs arising from the use of the products.

11. References

ISO 22174

Microbiology of food and animal feeding stuffs - Polymerase chain reaction (PCR) for the detection of food - borne pathogens - General requirements and definitions

ISO 20837

Microbiology of food and animal feeding stuffs - Polymerase chain reaction (PCR) for the detection of food - borne pathogens - Requirements for sample preparation and qualitative detection

ISO 20838

Microbiology of food and animal feeding stuffs - Polymerase chain reaction (PCR) for the detection of food - borne pathogens - Requirements for amplification and detection for qualitative methods

Khanna M., Fan J., Pehler-Harrington K., Waters C., Douglass P., Stallock J., Kehl S., Henrickson K.J.

“The Pneumoplex Assays, a Multiplex PCR-Enzyme Hybridization Assay That Allows Simultaneous Detection of Five Organisms, *M. pneumoniae*, *C.pneumoniae*, *Legionella pneumophila*, *Legionella micdadei*, and *B. pertussis*, and Its Real-Time Counterpart” *J. Clin. Microbio* *Legionella* (2005) 43, No. 2, p. 565–571.

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is a registered trademark of the ifp Institut für Produktqualität GmbH